XL t(12;21) ETV6/RUNX1 DF

Translocation/Dual Fusion Probe

Description

XL t(12;21) ETV6/RUNX1 DF is designed as a dual fusion probe. The orange labeled probe spans the breakpoint at 21q22 (RUNX1), the green labeled probe spans the breakpoint at 12p13 (ETV6).

XLt(12;21) ETV6/RUNX1 DF D-5115-100-OG is replacing XLt(12;21) D-5069-100-OG.

Clinical Details


Acute lymphoblastic leukemia (ALL) is a rapidly progressing cancer type characterized by the malignant transformation of lymphoid progenitor cells. It is the most common childhood cancer type and the second most leukemia in adults. Treatment of children usually results in good prognosis whereas the outcome for adults is less optimistic. Most patients show a transformation of precursors of the B-cell type, but also the Tcell phenotype is frequently observed. The most common aberration in pediatric Bcell ALL is t(12;21)(p13;q22) with an incidence of about 25% against <5% in adults. The result of this reciprocal translocation is an ETV6/RUNX1 fusion gene. Scientific data suggest that ETV6/RUNX1 is already established prenatally, but additional chromosomal aberrations are necessary for the development of ALL postnatally. The ETV6/RUNX1 fusion gene is transcriptionally active and is dysregulating a cascade of downstream genes. One study has shown, that all positive t(12;21) cases harbored the ETV6/RUNX1 fusion gene but not the reciprocal gene. This suggests, that ETV6/RUNX1 is involved in the manifestation of ALL, but not RUNX1/ETV6. Since t(12;21) is not detectable by conventional cytogenetic methods, FISH is one of the methods of choice.

Literature:

- **I** Romana et al (1995) Blood 85:3662-3670
- Leukemia 16:669-674
- Sun et al (2017) Oncotarget 8:35445-35459








XL t(12;21) ETV6/RUNX1 DF hybridized to lymphocytes. One normal and one aberrant interphase are shown. The expected normal signal pattern of XL t(12;21) ETV6/RUNX1 DF is two orange and two green signals, representing the two normal ETV6 and RUNX1 loci. The reciprocal translocation t(12;21)(p13;q22) is splitting the orange signal on chromsome 21 and the green signal on chromsome 12, resulting in two green/orange colocalization/fusion signals. The ETV6 loci on chromsome 12 is deleted. The remaining copy of RUNX1 is indicated by one orange signal. Deletion of the remaining ETV6 homolog is reported frequently.

Clinical Applications:

🛯 ALL



Related Products

Product	Size	Order No.
XL t(12;21), <i>REPLACED</i>	100µl	D-5069-100-0G
XL ETV6	100µl	D-5073-100-0G
XLRUNX1	100µl	D-5096-100-0G

MetaSystems Probes

MetaSystems Probes GmbH (Headquarters)

1. Industriestrasse 7 68804 Altlussheim, Germany tel +49 6205 2927 60 fax +49 6205 2927 29 info@metasystems-probes.com

MetaSystems Group, Inc.

70 Bridge Street Newton, MA 02458, USA tel +1 6179 2499 50 | fax +1 6179 2499 54 info@metasystems.org

MetaSystems S.r.l.

Via Gallarate 80 20151 Milano, Italy tel +39 0236 7587 51 fax +39 0245 3753 03 info@metasystems-italy.com

MetaSystems India Pvt., Ltd.

No. 1/1, 1st Floor, 1st Main Rd., 2nd cross Thimmaiah Garden, R T Nagar Bangalore Karnataka, 560 032, India tel +91 9535 7788 01 info@metasystems-india.com

MetaSystems Asia Co., Ltd.

Unit 108, 1/F, Bio-Informatics Centre No. 2 Science Park West Avenue Hong Kong Science Park Shatin, New Territories, Hong Kong tel +852 2587 8333 | fax +852 2587 8334 info@metasystems-asia.com

Document No. PFS-D5115-2017-08-01-S © 2017 by MetaSystems Probes

FACTSHEET



info@metasystems-probes.com www.metasystems-probes.com