

Description


The XL CDKN2A/9q22 locus-specific probe detects deletions in band 9p21. This probe is labeled in orange and hybridizes to CDKN2A (p16) as well as CDKN2B (p15). A green labeled probe hybridizes to a specific locus at 9q22 and functions as a reference probe.

Clinical Details

Acute lymphoblastic leukemia (ALL) is the most common childhood cancer type. T-cell acute lymphoblastic leukemia (T-ALL) is an aggressive and quickly progressing type of ALL affecting T-lymphocytes. Inactivation of the tumor suppressor genes CDKN2A/2B, located at chromosomal region 9p21, is a significant event in the development of T-ALL and other cancer types. Principally, loss of function might happen by deletion, methylation of promotor regions or mutations, whereas deletion of chromosomal region 9p21 seems to be the predominant mechanism. CDKN2A/2B deletions can be detected in about 60% of pediatric and about 50% of adult T-ALL cases. Most deletions are within the resolution of the FISH technique. Genetic alterations of the 9p21 locus result in loss of regulation of the cell cycle which is critical to cancer development.

Literature:

- Novara et al (2009) Human Genet 126:511-520
- Sulong et al (2009) Blood 113:100-107
- J Girardi et al (2017) Blood 129:1113-1123



Order No.:

D-5118-100-0G


XL CDKN2A/9q22 hybridized to lymphocytes. One normal interphase is shown. The expected signal pattern of XL CDKN2A/9q22, when hybridized to normal cells, is two orange and green seperated signals representing the two CDKN2A chromosomal regions and the respective controls. The locus specific control region labeled in green is generating a signal comparable to CDKN2A signals in size and intensity. Heterozygous deletions of CDKN2A are indicated by the loss of one orange signal, homozygous deletions by the loss of both orange signals per cell.

Clinical Applications:

📕 ALL







Related Products

Product	Size	Order No.
XL CDKN2A	100 µl	D-5053-100-0G

MetaSystems Probes

MetaSystems Probes GmbH (Headquarters)

1. Industriestrasse 7 68804 Altlussheim, Germany tel +49 6205 2927 60 fax +49 6205 2927 29 info@metasystems-probes.com

MetaSystems Group, Inc.

70 Bridge Street Newton, MA 02458, USA tel +1 6179 2499 50 | fax +1 6179 2499 54 info@metasystems.org

MetaSystems S.r.l.

Via Gallarate 80 20151 Milano, Italy tel +39 0236 7587 51 fax +39 0245 3753 03 info@metasystems-italy.com

MetaSystems India Pvt., Ltd.

No. 1/1, 1st Floor, 1st Main Rd., 2nd cross Thimmaiah Garden, R T Nagar Bangalore Karnataka, 560 032, India tel +91 9535 7788 01 info@metasystems-india.com

MetaSystems Asia Co., Ltd.

Unit 108, 1/F, Bio-Informatics Centre No. 2 Science Park West Avenue Hong Kong Science Park Shatin, New Territories, Hong Kong tel +852 2587 8333 | fax +852 2587 8334 info@metasystems-asia.com

Document No. PFS-D5118-2017-07-01-S © 2017 by MetaSystems Probes





info@metasystems-probes.com www.metasystems-probes.com