


Order No.: D-5123-100-OG

Description

XL 4q12 DC is a dual color probe which detects rearrangements and deletions at 4q12. The probe contains an orange probe hybridizing to the CHIC2 gene region. A green probe hybridizes proximal to FIP1L1, and another green probe hybridizes to PDGFRA and the region distal to the gene. XL 4q12 DC is a two color alternative to XL 4q12 D-5063-100-TC.

Clinical Details

The updated (2016) World Health Organization (WHO) classification of tumors of hematopoietic and lymphoid tissues indicates the category myeloid/lymphoid neoplasms with eosinophilia and rearrangement of PDGFRA, PDGFRB, FGFR1, or with PCM1-JAK2. Eosinophilia is a condition in which the number of eosinophilic granulocytes in peripheral blood or tissue is increased above the normal level. The hypereosinophilic syndrome (HES) is a rare hematologic disorder characterized by a marked and persistent increase of eosinophilic granulocytes in blood or tissue accompanied by organ damage. HES is associated with neoplastic (primary) or reactive (secondary) processes. Patients suffering from hematologic disorders with eosinophilia are typically diagnosed with chronic eosinophilic leukemia (CEL), myeloproliferative neoplasms (MPN), certain variants of acute myeloid leukemia, systemic mastocytosis and others. The most frequent aberration detected in CEL (10-20%) is the interstitial deletion of the CHIC2 gene with breakpoints in the FIP1L1 and PDGFRA genes. The deletion of a fragment of about 800kb is resulting in the FIP1L1-PDGFRA fusion gene, a constitutively activated tyrosine kinase, transforming hematopoietic cells. Other aberrations associated with CEL are translocations generating fusion genes involving PDGFRA. PDGFRA and PDGFRB fusions are sensitive to Imatinib, a proven treatment for BCR-ABL-positive chronic myeloid leukemia, whereas FGFR1 mutations are resistant.



XL 4q12 hybridized to lymphocytes, one normal interphase is shown. The expected normal signal pattern is two green-orange colocalization/fusion signals, representing two non rearranged FIP1L1-CHIC2-PDGFRA loci. An interstitial deletion of the CHIC2 gene is resulting in one green signal and one green-orange colocalization/fusion signal. PDGFRA translocations are indicated by two green-orange colocalization/fusion signal plus one seperated green signal.

Clinical Applications:


■ CML/MPN

Literature:

- Cools et al (2003) N Engl | Med 348:1201-1214
- Pardanani et al (2003) Blood 102:3093-3096
- Valent et al (2012) Expert Rev Hematol 5:157-176

FACTSHEET





Related Products

Product	Size	Order No.
XL4q12	100µl	D-5063-100-TC



MetaSystems Probes GmbH (Headquarters)

1. Industriestrasse 7 68804 Altlussheim, Germany tel +49 6205 2927 60 fax +49 6205 2927 29 info@metasystems-probes.com

MetaSystems Group, Inc.

70 Bridge Street Newton, MA 02458, USA tel +1 6179 2499 50 | fax +1 6179 2499 54 info@metasystems.org

MetaSystems S.r.l.

Via Gallarate 80 20151 Milano, Italy tel +39 0236 7587 51 fax +39 0245 3753 03 info@metasystems-italy.com

MetaSystems India Pvt., Ltd.

No. 1/1, 1st Floor, 1st Main Rd., 2nd cross Thimmaiah Garden, R T Nagar Bangalore Karnataka, 560 032, India tel +91 9535 7788 01 info@metasystems-india.com

MetaSystems Asia Co., Ltd.

Unit 108, 1/F, Bio-Informatics Centre No. 2 Science Park West Avenue Hong Kong Science Park Shatin, New Territories, Hong Kong tel +852 2587 8333 | fax +852 2587 8334 info@metasystems-asia.com

Document No. PFS-D5123-2017-08-01-S © 2017 by MetaSystems Probes

FACTSHEET

